




ANGLED BLADE PLATES & INSTRUMENTS



INSTRUMENTS SET

ABP-001 TRIPPLE DRILL GUIDE 130°



ABP-002 ROUTER



ABP-003 CHIESEL GUIDE



ABP-004 CEANG CHIESEL



ABP-005 INSERTER/EXTRACTER



ABP-006 PLATE IMPACTER



ABP-007 SLOTTED HAMMER



ABP-008 POSITIONING PLATE 90/50/40°

ABP-009 POSITIONING PLATE 80/70/30°

ABP-010 POSITIONING PLATE 100/60/20°



ABP-011 ANGLED PLATE GUIDE



AB-012 NEUTRAL & LOAD DRILL GUIDE




PLATES:

Cat. No. Holes Blade Length Shaft Length

CONDYLAR PLATES 95°

ABP-101	5	40MM	92MM
ABP-102	5	50MM	92MM
ABP-103	5	60MM	92MM
ABP-104	5	70MM	92MM
ABP-105	7	40MM	124MM
ABP-106	7	50MM	124MM
ABP-107	7	60MM	124MM
ABP-108	7	70MM	124MM
ABP-109	9	40MM	156MM
ABP-110	9	50MM	156MM
ABP-111	9	60MM	156MM
ABP-112	9	70MM	156MM
ABP-113	12	50MM	204MM
ABP-114	12	60MM	204MM
ABP-115	12	70MM	204MM
ABP-116	12	80MM	204MM
ABP-117	14	50MM	236MM
ABP-118	14	60MM	236MM
ABP-119	14	70MM	236MM
ABP-120	14	80MM	236MM
ABP-121	16	50MM	268MM
ABP-122	16	60MM	268MM
ABP-123	16	70MM	268MM
ABP-124	16	80MM	268MM
ABP-125	18	50MM	300MM
ABP-126	18	60MM	300MM
ABP-127	18	70MM	300MM
ABP-128	18	80MM	300MM



THICKNESS 5.7MM

WIDTH 16.2MM


HOLE SPACING 16.0MM

U PROFILE BLADE 6.6/16.2MM

IN TITANIUM: ABPT-101 TO ABPT-128

ANGLED BLADE PLATES 130°

ABP-301	4	50MM	60MM
ABP-302	4	60MM	60MM
ABP-303	4	70MM	60MM
ABP-304	4	80MM	60MM
ABP-305	4	85MM	60MM
ABP-306	4	90MM	60MM
ABP-307	4	95MM	60MM
ABP-308	4	100MM	60MM
ABP-309	4	105MM	60MM
ABP-310	4	110MM	60MM
ABP-311	6	50MM	104MM
ABP-312	6	60MM	104MM
ABP-313	6	70MM	104MM
ABP-314	6	80MM	104MM
ABP-315	6	90MM	104MM
ABP-316	9	50MM	152MM
ABP-317	9	60MM	152MM
ABP-318	9	70MM	152MM
ABP-319	9	80MM	152MM
ABP-320	9	90MM	152MM
ABP-321	12	70MM	200MM
ABP-322	12	80MM	200MM
ABP-323	12	90MM	200MM
ABP-324	14	70MM	232MM
ABP-325	14	80MM	232MM
ABP-326	14	90MM	232MM



THICKNESS 5.4MM
WIDTH 14.0MM
HOLE SPACING 16.0MM
T-PROFILE BLADE 5.0/11.7MM

IN TITANIUM: ABPT-301 TO ABPT-326

HIP PLATE 110°

ABP-401	4	65MM
ABP-402	4	75MM
ABP-403	4	85MM

THICKNESS 5.7MM
WIDTH 16.2MM
HOLE SPACING 16.0MM
U-PROFILE BLADE 6.6/16.2MM



IN TITANIUM: ABPT-401 TO ABPT-403

HIP PLATE 120°

ABP-501	4	65MM
ABP-502	4	75MM
ABP-503	4	85MM

THICKNESS 5.7MM
WIDTH 16.2MM
HOLE SPACING 16.0MM
U-PROFILE BLADE 6.6/16.2MM

IN TITANIUM: ABPT-501 TO ABPT-503

**HIP PLATE 130°**

ABP-601	4	65MM
ABP-602	4	70MM
ABP-603	4	75MM
ABP-603	4	80MM
ABP-601	4	85MM
ABP-602	4	90MM
ABP-603	4	95MM
ABP-601	4	100MM

THICKNESS 5.7MM
WIDTH 16.2MM
HOLE SPACING 16.0MM
U-PROFILE BLADE 6.6/16.2MM

IN TITANIUM: ABPT-601 TO ABPT-603



www.aonesurgico.com

Cat. No. Holes Displacement Shaft Length

CHILD HIP PLATE

ABP-701	3	8MM/80°	35MM
ABP-702	3	8MM/90°	35MM
ABP-703	3	8MM/100°	35MM

IN TITANIUM: ABPT-701 TO ABPT-703

THICKNESS 3.6MM
WIDTH 11.4MM
HOLE SPACING 16.0MM
T-PROFILE BLADE 4.5/11.3MM



ADULT HIP PLATE 90°

ABP-801	4	10MM	40MM
ABP-802	4	10MM	60MM

IN TITANIUM: ABPT-801 TO ABPT-802

THICKNESS 5.7MM
WIDTH 16.12MM
HOLE SPACING 16.0MM
U-PROFILE BLADE 6.5/16.2MM



www.aonesurgico.com