


**CANNULATED SCREW**  
**INSTRUMENTS & IMPLANTS**



CSL-001 SELF HOLDING SCREW FORCEP



CSL-002 PARRELEL WIRE GUIDE FOR GUIDEWIRES 2.0MM



CSL-003 CANNULATED SCREW DRIVER



CSL-004 HOLDING SLEEVE




CSL-005 DIRECT MEASURING GUIDE



CSL-006 CANNULATED DRILL BIT 4.5MM/230MM



CSL-007 CANNULATED TAP 7.0MM



CSL-009 THREADED GUIDE WIRE 2.0MM/230MM



**IMPLANTS:****7.0MM CANNULATED SCREWS WITH 16MM THREADS**

CSL-101	30MM
CSL-102	35MM
CSL-103	40MM
CSL-104	45MM
CSL-105	50MM
CSL-106	55MM
CSL-107	60MM
CSL-108	65MM
CSL-109	70MM
CSL-110	75MM
CSL-111	80MM
CSL-112	85MM
CSL-113	90MM
CSL-114	95MM
CSL-115	100MM
CSL-116	105MM
CSL-117	110MM
CSL-118	115MM
CSL-119	120MM



*HEAD DIAMETER 8.0MM*  
*HEXAGONAL SOCKET WIDTH 3.5MM*  
*SHAFT DIAMETER 4.5MM*  
*THREAD DIAMETER 7.0MM*  
*CORE DIAMETER 4.5MM*  
*THREAD LENGTH 16.0MM*  
*CANNULATION 2.1MM*  
*DRILL BIT 4.5MM*  
*TAP 7.0MM*

IN TITANIUM: CSLT-101 TO CSLT-119

**7.0MM CANNULATED SCREWS WITH 32MM THREADS**

CSL-301	30MM
CSL-302	35MM
CSL-303	40MM
CSL-304	45MM
CSL-305	50MM
CSL-306	55MM
CSL-307	60MM
CSL-308	65MM
CSL-309	70MM
CSL-310	75MM
CSL-311	80MM
CSL-312	85MM



*HEAD DIAMETER 8.0MM*  
*HEXAGONAL SOCKET WIDTH 3.5MM*  
*SHAFT DIAMETER 4.5MM*  
*THREAD DIAMETER 7.0MM*  
*CORE DIAMETER 4.5MM*  
*THREAD LENGTH 32.0MM*  
*CANNULATION 2.1MM*  
*DRILL BIT 4.5MM*



CSL-313	90MM	<i>TAP 7.0MM</i>
CSL-314	95MM	
CSL-315	100MM	
CSL-316	105MM	
CSL-317	110MM	
CSL-318	115MM	
CSL-319	120MM	

IN TITANIUM: CSLT-301 TO CSLT-319

**WASHER 13 MM DIAMETERS**

CSL-320



www.aonesurgico.com



## Cannulated Screw Instruments Set

IN STERILIZING CONTAINER

ITEM NUMBER: INS-O-CSS

CAT NO	DESCRIPTION	QUANTITY
CSS-001	SELF HOLDING SCREW FORCEPS	1
CSS-002	DOUBLE DRILL SLEEVE 2.7/1.2MM	1
CSS-003	SMALL CANNULATED SCREW DRIVER	1
CSS-004	HOLDING SLEEVE FOR CSS-O-003	1
CSS-005	DIRECT MEASURING DEVICE	1
CSS-006	CANNULATED DRILL BIT 2.7MM	1
CSS-007	CANNULATED TAP	1
CSS-008	STYLET 2.0MM	1
CSS-009	THREADED GUIDE WIRE 2.0MM	1
CSS-010	DRILL SLEEVE 3.5/2.7MM WITH STOP	1
CSS-011	COUNTERSINK CANNULATED	1

CSS-001 QUICK COUPLING TAP HANDLE

CSS-002 DOUBLE DRILL SLEEVE 2.7/1.2MM



CSS-003 SELF HOLDING SCREW FORCEPS



CSS-004 CANULATED SCREW DRIVER



CSS-005 HOLDING SLEEVE



CSS-006 DIRECT MEASURING DEVICE



CSS-007 DRILL SLEEVE WITH STOP 3.5/2.7MM



CSS-008 CANNULATED DRILL BIT 2.7MM/160MM



CSS-009 CANNULATED TAP 4.0MM



CSS-010 CANNULATED COUNTERSINK



CSS-011 THREADED GUIDE PIN 1.25/150MM



www.aonesurgico.com

**IMPLANTS:****4.0MM CANNULATED SCREWS SHORT THREADED**

CSS-101	10MM
CSS-102	12MM
CSS-103	14MM
CSS-104	16MM
CSS-105	18MM
CSS-106	20MM
CSS-107	22MM
CSS-108	24MM
CSS-109	26MM
CSS-110	28MM
CSS-111	30MM
CSS-112	32MM
CSS-113	34MM
CSS-114	36MM
CSS-115	38MM
CSS-116	40MM
CSS-117	42MM
CSS-118	44MM
CSS-119	46MM
CSS-120	48MM
CSS-121	50MM



*HEAD DIAMETER 5.0MM  
HEXAGONAL SOCKET WIDTH 2.5MM  
THREAD DIAMETER 4.0MM  
CORE DIAMETER 25MM  
SHAFT DIAMETER 2.7MM  
DRILL BIT 2.7MM  
TAP 4.0MM*

IN TITANIUM: CSST-101 TO CSST-121

**WASHER SMALL 7.0 MM DIAMETERS**

CSS-122

