




**MINI INSTRUMENTS & IMPLANTS**



WWW.AONESURGICO.COM

### JACOB CHUCK DRILL BITS

MIS-001	1.1MM	45MM LENGTH
MIS-002	1.5MM	70MM LENGTH
MIS-003	2.0MM	85MM LENGTH
MIS-004	2.7MM	85MM LENGTH



### QUICK COUPLING DRILL BITS

MIS-005	1.1MM	60MM LENGTH
MIS-006	1.5MM	85MM LENGTH
MIS-007	1.5MM	110MM LENGTH
MIS-008	2.0MM	100MM LENGTH
MIS-009	2.0MM	125MM LENGTH
MIS-010	2.7MM	100MM LENGTH
MIS-011	2.7MM	125MM LENGTH



### QUICK COUPLING TAP

MIS-012	1.5MM	50MM LENGTH
MIS-013	2.0MM	53MM LENGTH
MIS-014	2.0MM	100MM LENGTH
MIS-015	2.7MM	100MM LENGTH
MIS-016	2.7MM	125MM LENGTH



### QUICK COUPLING SCREW DRIVER SHAFT

MIS-017



### HEXAGONAL SCREW DRIVER WITH GROOVE

MIS-018



MIS-019 HOLDING SLEEVE



MIS-020 MINI SCREW DRIVER FOR 1.5MM-2.0MM SCREWS



MIS-021 COUNTERSINK



MIS-022 MINI QUICK COUPLING HANDLE



MIS-023 QUICK COUPLING HANDLE



MIS-024 DOUBLE DRILL SLEEVE 1.1MM/1.5MM



MIS-025 DOUBLE DRILL SLEEVE 1.5MM/2.0MM



MIS-026 NEUTRAL & LOAD DRILL SLEEVE  
FOR 1.5MM DRILL BITS



MIS-027 NEUTRAL & LOAD DRILL SLEEVE  
FOR 2.0MM DRILL BITS



MIS-028 SHORT DRILL GUIDE FOR 2.7MM PLATES



MIS-029 LONG DRILL GUIDE FOR 2.7MM PLATES



MIS-030 MINI DEPTH GAUGE FOR  
1.5MM-2.0MM SCREWS



MIS-031 SMALL DEPTH GAUGE FOR  
2.7MM-4.0MM SCREWS



www.aone.com

MIS-032 SHARP HOOK



MIS-033 SELF HOLDING SCEW FORCEPS



MIS-034 BENDING IRON FOR SMALL PLATES



MIS-035 SMALL RETRACTOR 6MM WIDE 160MM LENGTH  
SHORT NARROW TIP



MIS-037 SMALL RETRACTOR 15MM



MIS-038 PERIOSTEAL ELEVATOR 3MM WIDE



MIS-039 MINI PLATEHOLDING FORCEPS



MIS-040 SMALL REDUCTION FORCEPS POINTED  
130MM LENGTH



MIS-041 K-WIRES BENDING IRON 120MM LENGTH



MIS-042 MULTIPLE HOLE DRILL SLEEVE



MIS-043 WIRE BEIDING PLIER 155MM LENGTH



## 2.0MM MINI DYNAMIC COMPRESSION PLATES

CAT NO	HOLES	LENGTH
MCP-001	4	22MM
MCP-002	5	27MM
MCP-003	6	32MM
MCP-004	7	37MM
MCP-005	8	42MM



**THICKNESS 1.3MM**  
**WIDTH 5.1MM**  
**HOLE SPACING 5.0MM**

IN TITANIUM: MCPT-001 TO MCPT-005

## 2.0MM MINI L & R PLATES

MCP-101	LEFT
MCP-102	RIGHT
MCP-103	LEFT OBLIQUE
MCP-104	RIGHT OBLIQUE
MCP-105	T-SHAPED



**THICKNESS 1.0MM**  
**WIDTH 5.0MM**

IN TITANIUM MCPT-101 TO MCPT-105

## 1.5MM MINI STRAIGHT PLATE

MCP-201 20 100MM



IN TITANIUM: MCPT-201

## 2.0MM MINI ADAPTION PLATE

MCP-301 20 100MM



IN TITANIUM: MCPT-301

## 1.5MM MINI T-PLATES

MCP-401	3/3
MCP-402	3/5
MCP-403	4/4
MCP-404	4/6



IN TITANIUM: MCPT-401 TO MCP-404

## 2.0MM MINI T-PLATES

MCP-501	3/3
MCP-502	3/5
MCP-503	4/4
MCP-504	4/6



IN TITANIUM: MCPT-501 TO MCPT-504

## 2.7MM STRAIGHT RECONSTRUCTION PLATES

MCP-601	6	48MM
MCP-602	8	64MM
MCP-603	10	80MM
MCP-604	12	96MM
MCP-605	14	112MM
MCP-606	16	128MM



**THICKNESS 3.2MM**  
**WIDTH 8.0MM**  
**HOLE SPACING 8.0MM**

IN TITANIUM: MCPT-601 TO MCPT-606

## 2.7MM MINI DYNAMIC COMPRESSION PLATES

MCP-701	2	20MM
MCP-702	3	28MM
MCP-703	4	36MM
MCP-704	5	44MM
MCP-705	6	52MM
MCP-706	7	60MM
MCP-707	8	68MM
MCP-708	9	76MM
MCP-709	10	84MM
MCP-710	11	92MM
MCP-711	12	100MM



**THICKNESS 2.3MM**  
**WIDTH 8.1MM**  
**HOLE SPACING 8.0 & 12.0MM**

IN TITANIUM: MCPT-701 TO MCPT-711



### 1/4th TUBULAR PLATES

MCP-801	3	23MM
MCP-802	4	31MM
MCP-803	5	39MM
MCP-804	6	47MM
MCP-805	7	55MM
MCP-806	8	63MM



**THICKNESS 1.0MM**  
**WIDTH 7.0MM**  
**HOLE SPACING 8.0MM**

IN TITANIUM: MCPT-801 TO MCPT-808

### 2.7MM MINI L&T PLATES

MCP-901	LEFT
MCP-902	RIGHT
MCP-903	LEFT OBLIQUE
MCP-904	RIGHT OBLIQUE
MCP-905	T-SHAPED



**THICKNESS 1.2MM**  
**WIDTH 7.0MM**

IN TITANIUM: MCPT-901 TO MCPT-905

### K-WIRES

MCP-1001	0.8MM	70MM
MCP-1002	1.0MM	150MM
MCP-1003	1.2MM	150MM
MCP-1004	1.4MM	150MM
MCP-1005	1.5MM	150MM
MCP-1006	1.6MM	150MM
MCP-1007	1.8MM	150MM
MCP-1008	2.0MM	150MM
MCP-1009	2.5MM	150MM
MCP-1010	3.0MM	150MM

IN TITANIUM MCPT-1001 TO MCPT-1010



## SUTURE WIRES

MCP-1101	0.1MM	10 METERS
MCP-1102	0.2MM	10 METERS
MCP-1103	0.3MM	10 METERS
MCP-1104	0.4MM	10 METERS
MCP-1105	0.5MM	10 METERS



IN TITANIUM: MCPT-1101 TO MCP-1105

## 1.5MM CORTEX SCREWS

CAT NO	LENGTH
MSS-001	6MM
MSS-002	7MM
MSS-003	8MM
MSS-004	9MM
MSS-005	10MM
MSS-006	11MM
MSS-007	12MM
MSS-008	14MM
MSS-009	16MM
MSS-010	18MM
MSS-011	20MM



**HEAD DIAMETER 3.0MM**  
**HEXAGONAL SOCKET WIDTH 1.5MM**  
**THREAD DIAMETER 1.5MM**  
**CORE DIAMETER 1.0MM**  
**DRILL BIT FOR THREADED HOLES 1.1MM**  
**DRILL BIT FOR GLIDING HOLES 1.5MM**  
**TAP 1.5MM**

IN TITANIUM MSST-001 TO MSST-011

## 2.0MM CORTEX SCREWS

MSS-102	6MM
MSS-103	7MM
MSS-104	8MM
MSS-105	9MM
MSS-106	10MM
MSS-107	11MM
MSS-108	12MM
MSS-109	14MM
MSS-110	16MM
MSS-111	18MM
MSS-112	20MM
MSS-113	22MM



**HEAD DIAMETER 4.0MM**  
**HEXAGONAL SOCKET WIDTH 1.5MM**  
**THREAD DIAMETER 2.0MM**  
**CORE DIAMETER 1.3MM**  
**DRILL BIT FOR THREADED HOLES 1.5MM**  
**DRILL BIT FOR GLIDING HOLES 2.0MM**  
**TAP 2.0MM**

MSS-114	24MM
MSS-115	26MM
MSS-116	28MM
MSS-117	30MM
MSS-118	32MM
MSS-119	34MM
MSS-120	36MM
MSS-121	38MM
MSS-122	40MM

IN TITANIUM MSST-501 TO MSS-522

### 2.7MM CORTEX SCREWS

MSS-201	6MM
MSS-202	7MM
MSS-203	8MM
MSS-204	9MM
MSS-205	10MM
MSS-206	11MM
MSS-207	12MM
MSS-208	14MM
MSS-209	16MM
MSS-210	18MM
MSS-211	20MM
MSS-212	22MM
MSS-213	24MM
MSS-214	26MM
MSS-215	28MM
MSS-216	30MM
MSS-217	32MM
MSS-218	34MM
MSS-219	36MM
MSS-220	38MM
MSS-221	40MM



**HEAD DIAMETER 4.0MM**  
**HEXAGONAL SOCKET WIDTH 1.5MM**  
**THREAD DIAMETER 2.0MM**  
**CORE DIAMETER 1.3MM**  
**DRILL BIT FOR THREADED HOLES 1.5MM**  
**DRILL BIT FOR GLIDING HOLES 2.0MM**  
**TAP 2.0MM**

IN TITANIUM MSS-601 TO MSS-621

### SMALL WASHER 7.0MM DIAMETER

MSS-701

