

## **SMALL SCREW SET**



**3.5MM CORTEX SCREWS**

CAT NO	SIZE
SSS-001	10MM
SSS-002	12MM
SSS-003	14MM
SSS-004	16MM
SSS-005	18MM
SSS-006	20MM
SSS-007	22MM
SSS-008	24MM
SSS-009	26MM
SSS-010	28MM
SSS-011	30MM
SSS-012	32MM
SSS-013	34MM
SSS-014	36MM
SSS-015	38MM
SSS-016	40MM
SSS-017	42MM
SSS-018	44MM
SSS-019	46MM
SSS-020	48MM
SSS-021	50MM
SSS-022	55MM
SSS-023	60MM
SSS-024	65MM
SSS-025	70MM



**HEAD DIAMETER 6.0MM**  
**HEXAGONAL SOCKET WIDTH 2.5MM**  
**THREAD DIAMETER 3.5MM**  
**CORE DIAMETER 2.4MM**  
**DRILL BIT FOR GLIDING HOLES**  
**DRILL BIT FOR THREADED HOLES**

IN TITANIUM: SSS-001 TO SSS-025

**3.5MM CORETX SCREW SELF TAPPING**

SSS-201	10MM
SSS-202	12MM
SSS-203	14MM
SSS-204	16MM
SSS-205	18MM
SSS-206	20MM
SSS-207	22MM
SSS-208	24MM



SSS-209	26MM
SSS-210	28MM
SSS-211	30MM
SSS-212	32MM
SSS-213	34MM
SSS-214	36MM
SSS-215	38MM
SSS-216	40MM
SSS-217	42MM
SSS-218	44MM
SSS-219	46MM
SSS-220	48MM
SSS-221	50MM
SSS-222	55MM
SSS-223	60MM
SSS-224	65MM
SSS-225	70MM

**HEAD DIAMETER 6.0MM**  
**HEXAGONAL SOCKET WIDTH 2.5MM**  
**THREAD DIAMETER 3.5MM**  
**CORE DIAMETER 2.4MM**  
**DRILL BIT FOR GLIDING HOLES 3.55MM**  
**DRILL BIT FOR THREADED HOLES**

IN TITANIUM: SSST-201 TO SSS-225

#### 4.0MM CANCELLOUS SCREW SHORT THREADED

SSS-401	10MM
SSS-402	12MM
SSS-403	14MM
SSS-404	16MM
SSS-405	18MM
SSS-406	20MM
SSS-407	22MM
SSS-408	24MM
SSS-409	26MM
SSS-410	28MM
SSS-411	30MM
SSS-412	32MM
SSS-413	34MM
SSS-414	36MM
SSS-415	38MM
SSS-416	40MM
SSS-417	42MM
SSS-418	45MM



**HEAD DIAMETER 6.0MM**  
**HEXAGONAL SOCKET WIDTH 2.5MM**  
**THREAD DIAMETER 4.0MM**  
**CORE DIAMETER 1.9MM**  
**DRILL BIT FOR THREADED HOLES 2.0MM**  
**TAP 4.0MM**

SSS-419	50MM
SSS-420	55MM
SSS-421	60MM
SSS-422	65MM
SSS-423	70MM

IN TITANIUM: SSST-401 TO SSST-423

### 4.0MM CANCELLOUS SCREW FULLY THREADED

SSS-601	10MM
SSS-602	12MM
SSS-603	14MM
SSS-604	16MM
SSS-605	18MM
SSS-606	20MM
SSS-607	22MM
SSS-608	24MM
SSS-609	26MM
SSS-610	28MM
SSS-611	30MM
SSS-612	32MM
SSS-613	34MM
SSS-614	36MM
SSS-615	38MM
SSS-616	40MM
SSS-617	42MM
SSS-618	45MM
SSS-619	50MM
SSS-620	55MM
SSS-621	60MM
SSS-622	65MM
SSS-623	70MM



**HEAD DIAMETER 6.0MM**  
**HEXAGONAL SOCKET WIDTH 2.5MM**  
**THREAD DIAMETER 4.0MM**  
**CORE DIAMETER 1.9MM**  
**DRILL BIT FOR THREADED HOLES 2.0MM**  
**TAP 4.0MM**

IN TITANIUM: SSST-601 TO SSST-623

### WASHER SMALL 7.0MM SIZE


SSS-801

IN TITANIUM: SSST-801



## LOCKING SCREW Ø 2.7MM (HEAD LCP 2.4)

CAT NO	SIZE
SSS-001A	10MM
SSS-002A	12MM
SSS-003A	14MM
SSS-004A	16MM
SSS-005A	18MM
SSS-006A	20MM
SSS-007A	22MM
SSS-008A	24MM
SSS-009A	26MM
SSS-010A	28MM
SSS-011A	30MM
SSS-012A	32MM
SSS-013A	34MM
SSS-014A	36MM
SSS-015A	38MM
SSS-016A	40MM
SSS-017A	42MM
SSS-018A	44MM
SSS-019A	46MM
SSS-020A	48MM
SSS-021A	50MM
SSS-022A	55MM
SSS-023A	60MM
SSS-024A	65MM
SSS-025A	70MM




**THREAD DIAMETER 2.7MM**  
**DRILL BIT FOR THREADED HOLE 2.0MM**  
**CORE DIAMETER 2.1MM**  
**HEAD DIAMETER 3.5MM**

IN TITANIUM: SSST-001A TO SSST-025A

## LOCKING SCREW Ø 2.7MM (HEAD LCP3.5) SELF TAPPING

SSS-006C	20MM
SSS-007C	22MM
SSS-008C	24MM
SSS-009C	26MM
SSS-010C	28MM
SSS-011C	30MM
SSS-012C	32MM
SSS-013C	34MM




SSS-014C	36MM
SSS-015C	38MM
SSS-016C	40MM
SSS-017C	42MM
SSS-018C	44MM
SSS-019C	46MM
SSS-020C	48MM
SSS-021C	50MM
SSS-022C	55MM
SSS-023C	60MM

IN TITANIUM: SSST-006C TO SSST-023C

### LOCKING SCREW Ø 3.5MM SELF TAPPING

SSS-201A	10MM
SSS-202A	12MM
SSS-203A	14MM
SSS-204A	16MM
SSS-205A	18MM
SSS-206A	20MM
SSS-207A	22MM
SSS-208A	24MM
SSS-209A	26MM
SSS-210A	28MM
SSS-211A	30MM
SSS-212A	32MM
SSS-213A	34MM
SSS-214A	36MM
SSS-215A	38MM
SSS-216A	40MM
SSS-217A	42MM
SSS-218A	44MM
SSS-219A	46MM
SSS-220A	48MM
SSS-221A	50MM
SSS-222A	55MM
SSS-223A	60MM
SSS-224A	65MM
SSS-225A	70MM



**HEAD DIAMETER 5.0MM**  
**THREAD DIAMETER 3.5MM**  
**CORE DIAMETER 2.9MM**  
**DRILL BIT FOR THREADED HOLES 2.8MM**  
**TAP 3.5MM**

SSS-226A	75MM
SSS-227A	80MM
SSS-228A	85MM
SSS-229A	90MM

IN TITANIUM: SSST-201A TO SSST-229A

### LOCKING SCREW Ø 3.5MM SELF DRILLING

SSS-201B	10MM
SSS-202B	12MM
SSS-203B	14MM
SSS-204B	16MM
SSS-205B	18MM
SSS-206B	20MM
SSS-207B	22MM
SSS-208B	24MM
SSS-209B	26MM
SSS-210B	28MM
SSS-211B	30MM
SSS-212B	32MM
SSS-213B	34MM
SSS-214B	36MM
SSS-215B	38MM
SSS-216B	40MM
SSS-217B	42MM
SSS-218B	44MM
SSS-219B	46MM
SSS-220B	48MM
SSS-221B	50MM
SSS-222B	55MM
SSS-223B	60MM
SSS-224B	65MM
SSS-225B	70MM
SSS-226B	75MM
SSS-227B	80MM
SSS-228B	85MM
SSS-229B	90MM



**HEAD DIAMETER 5.0MM**  
**THREAD DIAMETER 3.5MM**  
**CORE DIAMETER 2.9MM**  
**DRILL BIT FOR THREADED HOLES 2.8MM**  
**TAP 3.5MM**

IN TITANIUM: SSST-201B TO SSST-229B