


CRANIO-MAXILLO-FACIAL SYSTEM (TITANIUM)

TITANIUM MINI PLATING SYSTEM



SELF FITTING SCREW

MANUFACTURED FROM HIGH QUALITY TITANIUM AND AVAILABLE IN CROSS HEAD SQUARE HEAD OR HEXAGONAL HEAD
SELF FITTING REALIZES ON EASIER AND FASTER

HANDLING AS THE USE OF THE SCREWDRIVER WITH A SCREW HOLDING DEVICE IS NOT NECESSARY

SELF FITTING REQUIRES A SCREW DRIVER BLADE WITH SELF FITTING FEATURES

SELF FITTING CAN AS WELL BE USED WITH COMMON SCREWDRIVER WITH OR WITH OUT SCREW HOLDING DEVICE, BUT THE SELF FITTING CHARACTER OF THE SCREW WILL THEN BE LOST.

ALL SCREWS IN THE SYSTEM ARE SELF TAPPING SCREWS



SELF DRILLING SCREW

MANUFACTURED FROM HIGH QUALITY TITANIUM ALLOY HAVE FOLLOWING CHARACTERISTICS
SELF DRILLING DOES NOT REQUIRE PREDRILLING AND THEREFORE REDUCES OPERATION TIME.

SELF DRILLING WITH **SELF FITTING** SCREW HEADS.

SELF DRILLING WITH A SHARP SCREW TIP FOR INSERTION.

SELF DRILLING IS SILVER COLOR FOR EASY DIFFERENTIATION.

ALL SCREWS IN SYSTEM ARE **SELF TAPPING**

FITTING DRILLS FOR EACH SYSTEM

DRILL	SCREW
1.0MMO	1.2MMO
	1.4MMO EMERGENCY
1.1MMO	1.5MMO
1.4MMO	1.7MMO
	1.9MM EMERGENCY
1.5MMO	2.0MMO
1.9MMO	2.3MMO
	2.5MMO EMERGENCY
2.0MMO	2.7MMO
	3.0MM EMERGENCY

TITANIUM SCREWS SELF TAPPING



STANDARD



CROSS HEAD
SELF FITTING



CROSS HEAD
SELF FITTING/DRILLING



HEXAGONAL
SELF FITTING



SQUARE HEAD
SELF FITTING

O×MM	CAT NO	CAT NO	CAT NO	CAT NO	CAT NO
2.0×4.0	CMF-2004A	CMF-2004B	CMF-2004C	CMF-2004D	CMF-2004E
2.0×5.0	CMF-2005A	CMF-2005B	CMF-2005C	CMF-2005D	CMF-2005E
2.0×6.0	CMF-2006A	CMF-2006B	CMF-2006C	CMF-2006D	CMF-2006E
2.0×7.0	CMF-2007A	CMF-2007B	CMF-2007C	CMF-2007D	CMF-2007E
2.0×8.0	CMF-2008A	CMF-2008B	CMF-2008C	CMF-2008D	CMF-2008E
2.0×9.0	CMF-2009A	CMF-2009B	CMF-2009C	CMF-2009D	CMF-2009E
2.0×10.0	CMF-2010A	CMF-2010B	CMF-2010C	CMF-2010D	CMF-2010E
2.0×11.0	CMF-2011A	CMF-2011B	CMF-2011C	CMF-2011D	CMF-2011E
2.0×12.0	CMF-2012A	CMF-2012B	CMF-2012C	CMF-2012D	CMF-2012E
2.0×13.0	CMF-2013A	CMF-2013B	CMF-2013C	CMF-2013D	CMF-2013E
2.0×14.0	CMF-2014A	CMF-2014B	CMF-2014C	CMF-2014D	CMF-2014E
2.0×15.0	CMF-2015A	CMF-2015B	CMF-2015C	CMF-2015D	CMF-2015E
2.0×16.0	CMF-2016A	CMF-2016B	CMF-2016C	CMF-2016D	CMF-2016E
2.0×17.0	CMF-2017A	CMF-2017B	CMF-2017C	CMF-2017D	CMF-2017E
2.0×18.0	CMF-2018A	CMF-2018B	CMF-2018C	CMF-2018D	CMF-2018E
2.0×19.0	CMF-2019A	CMF-2019B	CMF-2019C	CMF-2019D	CMF-2019E
2.0×20.0	CMF-2020A	CMF-2020B	CMF-2020C	CMF-2020D	CMF-2020E
2.0×21.0	CMF-2021A	CMF-2021B	CMF-2021C	CMF-2021D	CMF-2021E
2.0×22.0	CMF-2022A	CMF-2022B	CMF-2022C	CMF-2022D	CMF-2022E
2.0×23.0	CMF-2023A	CMF-2023B	CMF-2023C	CMF-2023D	CMF-2023E
2.0×24.0	CMF-2024A	CMF-2024B	CMF-2024C	CMF-2024D	CMF-2024E

www.aonesurgico.com

TITANIUM EMERGENCY SCREWS SELF TAPPING



STANDARD

CROSS HEAD
SELF FITTING



HEXAGONAL HEAD
SELF FITTING

HEXAGONAL HEAD
SELF FITTING

O×MM	CAT NO	CAT NO	CAT NO	CAT NO
2.3MM×4.0MM	CMF-2304A	CMF-2304B	CMF-2304C	CMF-2304D
2.3MM×5.0MM	CMF-2305A	CMF-2305B	CMF-2305C	CMF-2305D
2.3MM×6.0MM	CMF-2306A	CMF-2306B	CMF-2306C	CMF-2306D
2.3MM×7.0MM	CMF-2307A	CMF-2307B	CMF-2307C	CMF-2307D
2.3MM×8.0MM	CMF-2308A	CMF-2308B	CMF-2308C	CMF-2308D
2.3MM×9.0MM	CMF-2309A	CMF-2309B	CMF-2309C	CMF-2309D
2.3MM×10.0MM	CMF-2310A	CMF-2310B	CMF-2310C	CMF-2310D
2.3MM×11.0MM	CMF-2311A	CMF-2311B	CMF-2311C	CMF-2311D
2.3MM×12.0MM	CMF-2312A	CMF-2312B	CMF-2312C	CMF-2312D
2.3MM×13.0MM	CMF-2313A	CMF-2313B	CMF-2313C	CMF-2313D
2.3MM×14.0MM	CMF-2314A	CMF-2314B	CMF-2314C	CMF-2314D
2.3MM×15.0MM	CMF-2315A	CMF-2315B	CMF-2315C	CMF-2315D
2.3MM×16.0MM	CMF-2316A	CMF-2316B	CMF-2316C	CMF-2316D
2.3MM×17.0MM	CMF-2317A	CMF-2317B	CMF-2317C	CMF-2317D
2.3MM×18.0MM	CMF-2318A	CMF-2318B	CMF-2318C	CMF-2318D
2.3MM×19.0MM	CMF-2319A	CMF-2319B	CMF-2319C	CMF-2319D
2.3MM×20.0MM	CMF-2320A	CMF-2320B	CMF-2320C	CMF-2320D
2.3MM×21.0MM	CMF-2321A	CMF-2321B	CMF-2321C	CMF-2321D
2.3MM×22.0MM	CMF-2322A	CMF-2322B	CMF-2322C	CMF-2322D
2.3MM×23.0MM	CMF-2323A	CMF-2323B	CMF-2323C	CMF-2323D
2.3MM×24.0MM	CMF-2324A	CMF-2324B	CMF-2324C	CMF-2324D

MINI PLATES

STRAIGHT PLATES



2 HOLES
CMF-10102

4 HOLES
W/O BRIDGE
CMF-10104

4 HOLES
MEDIUM BRIDGE
CMF-10184

4 HOLES
WITH BRIDGE
CMF-10194

6 HOLES
W/O BRIDGE
CMF-10106

6 HOLES
WITH BRIDGE
CMF-10196

8 HOLES
W/O BRIDGE
CMF-10108

16HOLES
W/O BRIDGE
CMF-1016

L-PLATES



RIGHT W/O BRIDGE CMF-10104	LEFT W/O BRIDGE CMF-10224	RIGHT SHORT BRIDGE CMF-10254	LEFT SHORT BRIDGE CMF-10264	RIGHT WITH BRIDGE CMF-10274	LEFT WITH BRIDGE CMF-10294	RIGHT LONG BRIDGE CMF-10234	LEFT LONG BRIDGE CMF-10244
----------------------------------	---------------------------------	------------------------------------	-----------------------------------	-----------------------------------	----------------------------------	-----------------------------------	----------------------------------

Y-PLATES



W/O BRIDGE CMF-10705	WITH BRIDGE CMF-10795
-------------------------	--------------------------

DOUBLE Y-PLATES



W/O BRIDGE CMF-10606	WITH HOLE BRIDGE CMF-10607	WITH BRIDGE CMF-10696
-------------------------	-------------------------------	--------------------------



LONG BRIDGE CMF-10698	CMF-10698
--------------------------	-----------

NOSE PLATES



Y-PLATE CMF-10910

ORBITA PLATES




SEMICIRCULAR SHORT CMF-10806



SEMICIRCULAR LONG CMF-10808

T-PLATES



20 HOLES	40 HOLES	90°	90°	80° RIGHT	80° LEFT	90°	90°
W/O BRIDGE	W/O BRIDGE	W/O BRIDGE	WITH BRIDGE	WITH BRIDGE	WITH BRIDGE	W/O BRIDGE	WITH BRIDGE
CMF-10120	CMF-10140	CMF-10306	CMF-10375	CMF-10355	CMF-10365	CMF-10316	CMF-10314

L-PLATES



QUADRATIC/RECTANGULAR PLATES



INSTRUMENTATION

DRILL GUIDE WITH HANDLE
2.0MM SCREWS

WITH INSERT	SINGLE HANDLE
CMF-20318	CMF-20218



DEPTH GAUGE
19CM, 7"
CMF-20018




SPIRAL DRILL BIT
1.5MMO
7.0MM STOP×50MM CMF-15050
7.0MM STOP×70MM CMF-15070
22MM STOP ×50MM CMF-15950
22MM STOP ×70MM CMF-15970
22MM STOP × 105MM CMF-15105



TRANSBUCCALE DRILL GUIDE

COMPLETE FOR SCREW DRIVER W/O SCREW HOLDING DEVICE
3.0MMO CMF-30000
HANDLE WITH CMF-30010
SLEVE CMF-30030
CHEK-HOLDING CMF-30040
TROCAR CMF-30020
DRILL GUIDE

COMPLETE FOR SCREW DRIVER WITH SCREW HOLDING DEVICE
6.0MMO CMF-6000
HANDLE WITH SLEVE CMF-60010
CHEEK – HOLDING CMF-60020
TROCAR CMF-60030
DRILLGUIDE CMF-60040






PLATE AND WIRE CUTTING PLIER WITH TC
18CM, 7"
CMF-70418



PLATE BENDING PLIER WITH PIN
13.5CM, 5"
CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4"
CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313

PLATE HOLDING CLAMP
14CM, 5-1/2"
CMF-50013



PLATE HOLDING FORK
18CM, 7"
CMF-2218



PLATE HOLDING FORCEPS STRAIGHT
18CM, 7"
CMF-40018



PLATE HOLDING FORCEPS CURVED
18CM, 7"
CMF-45018





SCREW DRIVER	STANDARD	CROSS HEAD	HEXAGONAL HEAD	SQUARE HEAD
DESMANTABE COMPLETE	CMF-20118	CMF-20218	CMF-20318	CMF-20418
HANDLE	CMF-00000	CMF-00000	CMF-00000	CMF-00000
SCREW HOLDING DEVICE	CMF-20000	CMF-20000	CMF-20000	CMF-20000
BLADE	CMF-20100	CMF-20200	CMF-20300	CMF-20400




SCREWDRIVER DESMENTABLE STANDARD	CROSS HEAD	HEXAGONAL HEAD	SQUARE HEAD
COMPLETE SELF FITTING	CMF-202185	CMF-203185	CMF-204185
HANDLE	CMF-00000	CMF-00000	CMF-00000
BLADE	CMF-202005	CMF-203005	CMF-204005



SCREW DRIVER	STANDARD	CROSS HEAD	HEXAGONAL HEAD	SQUARE HEAD
COMPLETE	CMF-20118	CMF-20218	CMF-20318	CMF-20418
SCREW DRIVER		CMF-202185	CMF-203185	CMF-205185
SELF FITTING				


THE LOW PROFILE SYSTEM

TITANIUM SCREWS
SELF-TAPPING




OxMM	STANDARD	CROSS HEAD (SEFL FITTING)
1.5x4.0	CMF-15104	CMF-15204
1.5x5.0	CMF-15105	CMF-15205
1.5x6.0	CMF-15106	CMF-15206
1.5x7.0	CMF-15107	CMF-15207
1.5x8.0	CMF-15108	CMF-15208
1.5x9.0	CMF-15109	CMF-15209
1.5x10.0	CMF-15110	CMF-15210
1.5x11.0	CMF-15111	CMF-15211
1.5x12.0	CMF-15112	CMF-15212
1.5x13.0	CMF-15113	CMF-15213
1.5x14.0	CMF-15114	CMF-15214
1.5x15.0	CMF-15115	CMF-15215
2.0x4.0	CMF-20704	CMF-20204
2.0x5.0	CMF-20705	CMF-20205
2.0x6.0	CMF-20706	CMF-20206
2.0x7.0	CMF-20707	CMF-20207
2.0x8.0	CMF-20708	CMF-20208
2.0x9.0	CMF-20709	CMF-20209
2.0x10.0	CMF-20710	CMF-20210
2.0x11.0	CMF-20711	CMF-20211
2.0x12.0	CMF-20712	CMF-20212
2.0x13.0	CMF-20713	CMF-20213
2.0x14.0	CMF-20714	CMF-20214
2.0x15.0	CMF-20715	CMF-20215
2.0x16.0	CMF-20716	CMF-20216
2.0x17.0	CMF-20717	CMF-20218
2.0x18.0	CMF-20718	CMF-20218
2.0x19.0	CMF-20719	CMF-20219
2.0x20.0		CMF20220
2.0x21.0		CMF-20221
2.0x22.0		CMF-20222
2.0x23.0		CMF-20223
2.0x24.0		CMF-20224

LOW PROFILE SYSTEM



HEIGHT 1.1MM THICKNESS 0.6MM

MINI PLATING SYSTEM



HEIGHT 1.3MM THICKNESS 1.0MM

LOW PROFILE PLATES

STRAIGHT PLATES



2 HOLES CMF-06102	4 HOLES W/O BRIDGE CMF-06104	4 HOLES MEDIUM BRIDGE CM F-06184	4 HOLES WITH BRIDGE CMF-06194	6 HOLES W/O BRIDGE CMF-06106	6 HOLES WITH BRIDGE CMF-06196	8 HOLES W/O BRIDGE CMF-06108
----------------------	------------------------------------	--	-------------------------------------	------------------------------------	-------------------------------------	------------------------------------

L-PLATES



RIGHT W/O BRIDGE CMF-06204	LEFT W/O BRIDGE CMF-06224	RIGHT SHORT BRIDGE CMF-06254	LEFT SHORT BRIDGE CMF-06264	RIGHT WITH BRIDGE CMF-06274	LEFT WITH BRIDGE CMF-06294	RIGHT LONG BRIDGE CMF-06275	LEFT LONG BRIDGE CMF-06295
----------------------------------	---------------------------------	------------------------------------	-----------------------------------	-----------------------------------	----------------------------------	-----------------------------------	----------------------------------

T-PLATES

DOUBLE T-PLATES



ORBITA PLATES



LOW PROFILE PLATES



Y-PLATES

DOUBLE Y-PLATES



NOSE PLATES



H-PLATE
CMF-06912



Y-PLATE
CMF-06910



BURR HOLE COVER
O 7MM CMF-06990
O 10MM CMF-06995

O 14MM CMF-06997

QUADRATIC AND RECTANGULAR PLATES



RECTANGULAR
QUADRATIC
CMF-06404
CMF-06504




CMF-06406
CMF-06506



CMF-06408
CMF-06508




CMF-06416
CMF-06416




CMF-06424
CMF-06524



14x2HOLES
CMF-06528



9x14 HOLES
CMF-06599



35x12 HOLES
CMF-06550

INSTRUMENTATION

DRILL GUIDE WITH HANDLE
2.0MM SCREWS

WITH INSERT SINGLE HANDLE
CMF-20318 CMF-20218



DEPTH GAUGE
19CM, 7"
CMF-20018



SPIRAL DRILL BIT
1.5MMO

7.0MM STOP×50MM CMF-15050
7.0MM STOP×70MM CMF-15070
22MM STOP ×50MM CMF-15950
22MM STOP ×70MM CMF-15970z
22MM STOP × 105MM CMF-15105



SPIRAL DRILL BIT
1.1MM O

5MM STOP×45MM CMF-11005
7MM STOP×45MM CMF-11007
9MM STOP×45MM CMF-11009
15MM STOP×45MM CMF-11015



TRANSBUCCALE DRILL GUIDE

COMPLETE FOR SCREW DRIVER W/O SCREW HOLDING DEVICE
3.0MMO CMF-30000
HANDLE WITH SLEVE CMF-30010
CHEK-HOLDING CMF-30020
TROCAR CMF-30030
DRILL GUIDE CMF-30040

COMPLETE FOR SCREW DRIVER WITH SCREW HOLDING DEVICE
6.0MMO CMF-60000
HANDLE WITH SLEVE CMF-60010
CHEEK – HOLDING CMF-60020
TROCAR CMF-60030
DRILLGUIDE CMF-60040






PLATE AND WIRE CUTTING PLIER WITH TC
18CM, 7" CMF-70418



PLATE BENDING PLIER WITH PIN
13.5CM, 5" CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4" CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313

PLATE HOLDING CLAMP
14CM, 5-1/2"
CMF-50013



PLATE HOLDING FORK
18CM, 7"
CMF-2218



PLATE HOLDING FORCEPS STRAIGHT
18CM, 7"
CMF-40018



PLATE HOLDING FORCEPS CURVED
18CM, 7"
CMF-45018





SCREW DRIVER	STANDARD	CROSS HEAD	HEXAGONAL HEAD	SQUARE HEAD
DESMANTABLE COMPLETE	CMF-20118	CMF-20218	CMF-20318	CMF-20418
HANDLE	CMF- 00000	CMF -00000	CMF- 00000	CMF- 00000
SCREW HOLDING DEVICE	CMF-20000	CMF- 20000	CMF-20000	CMF- 20000
BLADE	CMF-20100	CMF-20200	CMF-20300	CMF-20400



SCREWDRIVER DESMENTABLE	STANDARD	CROSS HEAD	HEXAGONAL HEAD	SQUARE HEAD
COMPLETE SELF FITTING		CMF-202185	CMF-203185	CMF-204185
HANDLE		CMF-00000	CMF- 00000	CMF- 00000
BLADE		CMF-202005	CMF-203005	CMF-204005



SCREW DRIVER	STANDARD	CROSS HEAD	HEXAGONAL HEAD	SQUARE HEAD
COMPLETE	CMF-20118	CMF-20218	CMF-20318	CMF-20418
SCREW DRIVER		CMF-202185	CMF-203185	CMF-205185
SELF FITTING				

CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

MEDIUM PLATING SYSTEM

TITANIUM SCREWS SELF TAPPING



OxMM	STANDARD CAT NO	CROSS HEAD SELF FITTING CAT NO	CROSS HEAD SELF FITTING/DRILLING CAT NO	SQUARE HEAD SELF FITTING CAT NO
1.5x3.0	CMF-15103	CMF-15203	CMF-15303	CMF-15403
1.5x4.0	CMF-15104	CMF-15204	CMF-15304	CMF-15404
1.5x5.0	CMF-15105	CMF-15205	CMF-15305	CMF-15405
1.5x6.0	CMF-15106	CMF-15206	CMF-15306	CMF-15406
1.5x7.0	CMF-15107	CMF-15207	CMF-15307	CMF-15407
1.5x8.0	CMF-15108	CMF-15208	CMF-15308	CMF-15408
1.5x9.0	CMF-15109	CMF-15209	CMF-15309	CMF-15409
1.5x10.0	CMF-15110	CMF-15210	CMF-15310	CMF-15410
1.5x11.0	CMF-15111	CMF-15211	CMF-15311	CMF-15411
1.5x12.0	CMF-15112	CMF-15212		CMF-15412
1.5x13.0	CMF-15113	CMF-15213		CMF-15413
1.5x14.0	CMF-15114	CMF-15214		CMF-15414
1.5x15.0	CMF-15115	CMF-15215		CMF-15415
1.5x16.0	CMF-15116	CMF-15216		CMA-15416
1.5x17.0	CMF-15117	CMF-15217		CMF-15417
1.5x18.0	CMF-15118	CMF-15218		CMF-15417
1.5x19.0	CMF-15119	CMF-15219		CMF-15419
1.5x20.0	CMF-15120	CMF-15220		CMF-15420
1.5x21.0		CMF-15221		
1.5x22.0		CMF-15222		
1.5x23.0		CMF-15223		
1.5x24.0		CMF-15224		

www.aonesurgico.com

TITANIUM EMERGENCY SCREWS
SELF TAPPING



OxMM	STANDARD CAT NO	CROSS HEAD SELF FITTING CAT NO	CROSS HEAD SELF FITTING/DRILLING CAT NO	SQUARE HEAD SELF FITTING CAT NO
1.7x3.0	CMF-17103	CMF-17203	CMF-17303	CMF-17403
1.7x4.0	CMF-17104	CMF-17204	CMF-17304	CMF-17404
1.7x5.0	CMF-17105	CMF-17205	CMF-17305	CMF-17405
1.7x6.0	CMF-17106	CMF-17206	CMF-17306	CMF-17406
1.7x7.0	CMF-17107	CMF-17207	CMF-17307	CMF-17407
1.7x8.0	CMF-17108	CMF-17208	CMF-17308	CMF-17408
1.7x9.0	CMF-17109	CMF-17209	CMF-17309	CMF-17409
1.7x10.0	CMF-17110	CMF-17210	CMF-17310	CMF-17410
1.7x11.0	CMF-17111	CMF-17211	CMF-17311	CMF-17411
1.7x12.0	CMF-17112	CMF-17212		CMF-17412
1.7x13.0	CMF-17113	CMF-17213		CMF-17413
1.7x14.0	CMF-17114	CMF-17214		CMF-17414
1.7x15.0	CMF-17115	CMF-17215		CMF-17415
1.7x16.0	CMF-17116	CMF-17216		CMA-17416
1.7x17.0	CMF-17117	CMF-17217		CMF-17417
1.7x18.0	CMF-17118	CMF-17218		CMF-17417
1.7x19.0	CMF-17119	CMF-17219		CMF-17419
1.7x20.0	CMF-17120	CMF-17220		CMF-17420
1.7x21.0	CMF-17121	CMF-17221		
1.7x22.0	CMF-17122	CMF-17222		
1.7x23.0	CMF-17123	CMF-17223		
1.7x24.0	CMF-17124	CMF-17224		
1.7x25.0	CMF-17125	CMF-17225		


TITANIUM EMERGENCY SCREWS
SELF TAPPING



	STANDARD	CROSS HEAD SELF FITING	SQUARE HEAD SELF FITING
O×MM	CAT NO	CAT NO	CAT NO
1.9×4.0	CMF-19104	CMF-19204	CMF-19304
1.9×5.0	CMF-19105	CMF-19205	CMF-19305
1.9×6.0	CMF-19106	CMF-19206	CMF-19306
1.9×7.0	CMF-19107	CMF-19207	CMF-19307
1.9×8.0	CMF-19108	CMF-19208	CMF-19308
1.9×9.0	CMF-19109	CMF-19209	CMF-19309
1.9×10.0	CMF-19110	CMF-19210	CMF-19310
1.9×11.0	CMF-19111	CMF-19211	CMF-19311
1.9×12.0	CMF-19112	CMF-19212	CMF-19312
1.9×13.0	CMF-19113	CMF-19213	CMF-19313
1.9×14.0	CMF-19114	CMF-19214	CMF-19314
1.9×15.0	CMF-19115	CMF-19215	CMF-19315
1.9×16.0	CMF-19116	CMF-19216	CMF-19316
1.9×17.0	CMF-19117	CMF-19217	CMF-19317
1.9×18.0	CMF-19118	CMF-19218	CMF-19318
1.9×19.0	CMF-19119	CMF-19219	CMF-19319
1.9×20.0	CMF-19120	CMF-19220	CMF-19320

www.aonesurgico.com

MEDIUM PLATES



STRAIGHT PLATES



2 HOLES
CMF-06102

4 HOLES
CMF-06104

6 HOLES
CMF-06106

8 HOLES
CMF-06108

12 HOLES
CMF-06112

16 HOLES
CMF-06116



4 HOLES
SHORT BRIDGE
CMF-06174

4 HOLES
MEDIUM BRIDGE
CMF-06184

4 HOLES
LONG BRIDGE
CMF-06194

6 HOLES
WITH BRIDGE
CMF-06196

24 HOLES
W/O BRIDGE
CMF-06124

30 HOLES
WITH BRIDGE
CMF-06130

L-PLATES



RIGHT 90° W/O BRIDGE CMF-06204	LEFT 90° W/O BRIDGE CMF-06224	RIGHT 90° SHORT BRIDGE CMF-06254	LEFT 90° SHORT BRIDGE CMF-06264	RIGHT 90° MED BRIDGE CMF-06274	LEFT 90° MEDI BRIDGE CMF-06294	RIGHT 90° LARGE BRIDGE CMF-06234	LEFT 90° LARGE BRIDGE CMF-06244
--------------------------------------	-------------------------------------	--	---------------------------------------	--------------------------------------	--------------------------------------	--	---------------------------------------



RIGHT EX LONG BRIDGE CMF-06275	LEFT EX LONG BRIDGE CMF-06295	RIGHT 112° W/O BRIDGE CMF-06210	LEFT 112° W/O BRIDGE CMF-06211	RIGHT 112° MEDIUM BRIDGE CMF-06214	LEFT 112° MEDIUM BRIDGE CMF-06215
--------------------------------------	-------------------------------------	---------------------------------------	--------------------------------------	--	---



RIGHT 112° LONG BRIDGE CMF-06216	LEFT 112° LONG BRIDGE CMF-06217	RIGHT 112° EXT LONG BRIDGE CMF-06218	LEFT 112° EXT LONG BRIDGE CMF-06219
--	---------------------------------------	--	---

T-PLATES

DOUBLE- T- PLATE

Y-PLATES



4 HOLES SHORT BRIDGE CMF-06304	4 HOLES MEDIUM BRIDGE CMF-06305	6 HOLES W/O BRIDGE CMF-06306	5 HOLES WITH BRIDGE CMF-0606307	6 HOLES W/O BRIDGE CMF-06316	5 HOLES W/O BRIDGE CMF-06705	5 HOLES MEDIUM BRIDGE CMF-06785	5 HOLES LONG BRIDGE CMF-06795
--------------------------------------	---------------------------------------	------------------------------------	---------------------------------------	------------------------------------	------------------------------------	---------------------------------------	-------------------------------------

DOUBLE- Y- PLATES



6 HOLES
W/O BRIDGE
CMF-06606

6 HOLES
WITH BRIDGE
CMF-06696

7 HOLES
W/O BRIDGE
CMF-06607

ORBITA PLATES



4 HOLES
W/O BRIDGE
CMF-06804

4 HOLES
SHORT BRIDGE
CMF-06814



10 HOLES W/O BRIDGE
CMF-06810



6 HOLES
MEDIUM BRIDGE
CMF-06816

8 HOLES
W/O BRIDGE
CMF-06808

BURR HOLE COVER



7MM O
CMF-06997



10MMO
CMF-06995



14MM O
CMF-06990

ZYGOMA PLATE



CMF-06518

QUADRATIC AND RECTANGULAR GRILLE-PLATES




2x2 HOLES
CMF-06504




3x2 HOLES
CMF-O-06506



4x2 HOLES
CMF-O-06508



4x4 HOLES
CMF-O-06516



6x4 HOLES
CMF-O-06524



2x2 HOLES
CMF-O-06404



3x2HOLES
CMF-06406



4x2 HOLES
CMF-06408



4x4 HOLES
CMF-06416



6x4 HOLES
CMF-06424

SPIRAL DRILL BIT

- 1.4MM Ø
- 5.0MM STOP×45MM CMF-14005
- 7.0MM STOP×45MM CMF-14007
- 9.0MM STOP 45M CMF-14009
- 15MM STOP ×45MM CMF-14015



SPIRAL DRILL BIT

- 1.1MM Ø
- 5MM STOP×45MM CMF-11005
- 7MM STOP×45MM CMF-11007
- 9MM STOP×45MM CMF-11009
- 15MM STOP×45MM CMF-11015



PLATE HOLDING FORCEPS STRAIGHT

18CM, 7" CMF-40018



PLATE HOLDING FORCEPS CURVED

18CM, 7" CMF-45018






PLATE AND WIRE CUTTING PLIER WITH TC
18CM, 7"
CMF-70418



PLATE BENDING PLIER WITH PIN
13.5CM, 5"
CMF-70213



PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4"
CMF-70112



FLAT PLIER
13.5CM, 5"
CMF-70313

PLATE HOLDING CLAMP
14CM, 5-1/2"
CMF-50013



SCREW DRIVER
DESAMTABE COMPLETE
HANDLE
SCREW HOLDING DEVICE
BLADE

STANDARD
CMF-20117
CMF-00000
CMF-20000
CMF-20100

CROSS HEAD
CMF-20217
CMF- 00000
CMF- 20000
CMF-20200

SQUARE HEAD
CMF-20417
CMF- 00000
CMF- 20000
CMF-20400



SCREWDRIVER DESMENTABLE
COMPLETE SELF FITTING
HANDLE
BLADE

STANDARD

CROSS HEAD
CMF-20217
CMF-00000
CMF-20200

SQUARE HEAD
CMF-20417
CMF- 00000
CMF-20400



SCREW DRIVER
COMPLETE
SCREW DRIVER
SELF FITTING

STANDARD
CMF-20117

CROSS HEAD
CMF-20217
CMF-20217

SQUARE HEAD
CMF-20417
CMF-20517

CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

MICRO PLATTING SYSTEM
TITANIUM SCREWS
SELF-TAPPING



STANDARD

CROSS HEAD
SELF FITTING

O×MM	CAT NO	CAT NO
1.2×2.0	CMF-12102	CMF-12202
1.2×3.0	CMF-12103	CMF-12203
1.2×4.0	CMF-12104	CMF-12204
1.2×5.0	CMF-12105	CMF-12205
1.2×6.0	CMF-12106	CMF-12206
1.2×7.0	CMF-12107	CMF-12207
1.2×8.0	CMF-12108	CMF-12208
1.2×9.0	CMF-12109	CMF-12209
1.2×10.0	CMF-12110	CMF-12210
1.2×11.0	CMF-12111	CMF-12211
1.2×12.0	CMF-12112	CMF-12212
1.2×13.0	CMF-12113	CMF-12213
1.2×14.0	CMF-12114	CMF-12214
1.2×15.0	CMF-12115	CMF-12215

TITANIUM EMERGENCY SCREWS
SELF-TAPPING




STANDARD

CROSS HEAD
SELF FITTING

O×MM	CAT NO	CAT NO
1.4×3.0	CMF-O-14103	CMF-O-14203
1.4×4.0	CMF-O-14104	CMF-O-14204
1.4×5.0	CMF-O-14105	CMF-O-14205
1.4×6.0	CMF-O-14106	CMF-O-14206
1.4×7.0	CMF-O-14107	CMF-O-14207
1.4×8.0	CMF-O-14108	CMF-O-14208
1.4×9.0	CMF-O-14109	CMF-O-14209
1.4×10.0	CMF-O-14110	CMF-O-14210
1.4×11.0	CMF-O-14111	CMF-O-14211
1.4×12.0	CMF-O-14112	CMF-O-14212
1.4×13.0	CMF-O-14113	CMF-O-14213
1.4×14.0	CMF-O-14114	CMF-O-14214
1.4×15.0	CMF-O-14115	CMF-O-14215

MICRO PLATES



STRAIGHT PLATES



2 HOLES
CMF-06102

4 HOLES
CMF-06104

6 HOLES
CMF-06106

8 HOLES
CMF-06108

12 HOLES
CMF-06112

16 HOLES
CMF-06116

24 HOLES
CMF-06124

L-PLATES



T-PLATES



5 HOLES
RIGHT SHORT L
CMF-06205

5 HOLES
LEFT SHORT
CMF-06245

8 HOLES
RIGHT LONG
CMF-0606208

8 HOLES
LEFT 8 HOLES
CMF-06248

90° 5 HOLES
SHORT
CMF-06305

90° 5 HOLES
LONG
CMF-06307

Y-PLATES



DOUBLE- Y-PLATES



H-PLATES



6 HOLES
CMF-06706


SHORT 6 HOLES
CMF-06606

LONG 7 HOLES
CMF-06607

7 HOLES
CMF-06807

9 HOLES
CMF-06809

14 HOLES
CMF-06814



QUADRATIC AND RECTANGULAR GRILLE-PLATES




4 HOLES
CMF-06504




6 HOLES
CMF-06506



8 HOLES
CMF-06508



16 HOLES
CMF-06516



24 HOLES
CMF-06524



4 HOLES
CMF-06404



6 HOLES
CMF-06406



8 HOLES
CMF-06408



16 HOLES
CMF-06416



24 HOLES
CMF-06424

BURR HOLE COVER



7MM O
CMF-06997



10MM O
CMF-06995



14MM
CMF-06990

ORBITA PLATE



SEMICIRCULAR, 8 HOLES
CMF-06901

TITANIUM MICRO MESH



0.1MM

CMF-1990	90×20MM
CMF-1991	120×60MM
CMF-1992	60×60MM

0.2MM

CMF-02990	90×120MM
CMF-02991	120×60MM
CMF-02992	60×60MM

SMITH WIRE CUTTING SCISSORS

TC. 16CM, 6"

STRAIGHT	CURVED
CMF-84416	CMF-84516



WIRE CUTTING SCISSOR TC

12CM

CMF-83512



INSTRUMENTATION

SPIRAL DRILL BIT

- 1.0MMO
- 4.0MM STOP×46MM CMF-15906
- 8.0MM STOP×50MM CMF-15907
- 12.0MM STOP 54M CMF-15908



DEPTH GAUGE

- 19CM, 7"
- CMF-20018



PLATE HOLDING FORCEPS STRAIGHT

- 18CM, 7" CMF-40018



PLATE HOLDING FORCEPS CURVED

- 18CM, 7" CMF-45018



PLATE AND WIRE CUTTING PLIER W

- 12.5CM, 4-3/4" CMF-70412



PLATE BENDING PLIER WITH PIN
13.5CM, 5" CMF-70213




PLATE BENDING PLIER WITH, ADEREA
12.5CM, 4-3/4" CMF-70112



FLAT PLIER
13.5CM, 5" CMF-70313



PLATE HOLDING CLAMP
14CM, 5-1/2" CMF-50013



DESMANTABLEABLE SCREW DRIVER
COMPLETE
STANDARD QCROSS HEAD
CMF-20117 CMF-20217
HANDLE
CMF-00000 CMF-00000
TENTION PLIER
CMF-20000 CMF-20000
BLADE
CMF-20100 CMF-20200



SCREW DRIVER DESMENTABLE
COMPLETE SELF FITTING
STANDARD CROSS HEAD
CMF-20117 CMF-20217
HANDLE CRF-00000
BLADE CMF-20200



SCREW DRIVER FIXED COMPLETE
STANDARD CROSS HEAD
CMF-20117 CMF-20217
SELF FITTING COMPLETE
SCREW DRIVER CMF-20217



CRANIO-MAXILLO-FACIAL TITANIUM SYSTEM

THE MANDIBULAR RECONSTRUCTION SYSTEM

TITANIUM SCREWS
SELF TAPPING

TITANIUM EMERGENCY
SCREW SELF TAPPING

TITANIUM FASTENING SCREWS




CROSS HEAD O×MM	CAT NO
2.7×8	CMF-27298
2.7×10	CMF-27210
2.7×12	CMF-27212
2.7×14	CMF-27214
2.7×16	CMF-27216
2.7×18	CMF-27218
2.7×20	CMF-27220
2.7×22	CMF-27222

CROSS HEAD O×MM	CAT NO
3.0×10	CMF-30210
3.0×12	CMF-30212
3.0×14	CMF-30214

CROSS HEAD CAT NO
CMF-28307

RECONSTRUCTION PLATES



8 HOLES
CMF-22008
CMF-22006



11 HOLES
CMF-22011
CMF-22011A



14 HOLES
CMF-22014
CMF-22014A



17 HOLES
CMF-22017
CMF-22017A

RECONSTRUCTION PLATE
TIN BENDING TEMPLATE



4+8 HOLES LEFT
CMF-22312
CMF-22312



4+8 HOLES RIGHT
CMF-22212
CMF-22212

RECONSTRUCTION PLATE
TIN BENDING TEMPLATE



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+11 HOLES LEFT
CMF-22315
CMF-22215

4+11 HOLES RIGHT
CMF-22215
CMF-22215



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE


4+14 HOLES LEFT
CMF-22318
CMF-22218

4+14 HOLES RIGHT
CMF-22218
CMF-22218



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+17+4 HOLES RIGHT
CMF-22425
CMF-22425



RECONSTRUCTION PLATE
TIN BENDING TEMPLATE

4+17+4 HOLES LEFT
CMF-22427
CMF-22427

PROSTHESES



TEMPORARY CONDYLAR PROSTHESES
RIGHT SHORT. CMF-25153




TEMPORARY CONDYLAR PROSTHESES
LEFT SHORT CMF-22150




TEMPORARY CONDYLAR
PROSTHESES LEFT 5 HOLES
CMF-25152



TEMPORARY CONDYLAR
PROSTHESES RIGHT 5 HOLES
CMF-25155



BENDING TEMPLATES CMF-22000



INSTRUMENTATION

PLATE BENDER 15 CM LENGTHS
CMF-22150




PLATE HOLDING FORCEPS LEFT
15.5CM CMF-25154



PLATE BENDING FORCEPS RIGHT
15.5 CM CMF-25155



PLATE BENDING PLIER
LENGTH 16.5MM FOR
2.2&2.5MM PLATES
CMF-25160





SCREWDRIVER DESMENTABLE
COMPLETE SELF FITTING
CROSS HEAD CMF-27218
HANDLE CMF-00000
BLADE CMF-27200



SCREWDRIVER DESMENTABLE
COMPLETE CMF-27218
HANDLE CMF-00000
BLADE CMF-27200
SCREW HOLDER CMF-27000A



TAP. 2.7MM SCREWS
CMF-27000



DRILL GUIDE (MANDIBULAR SYSTEM)
CMF-28100



SPIRAL DRILL O 2.0MM, 35 MM TOP
103MM CMF-20103

www.aonesurgico.com

DEPTH GAUGE (UP TO 45 MM SCREWS)
CMF-27018



PLATE CUTTING PLIERS 57.5CM LENGTH
CMF-47000



TRANSBUCCALE DRILL GUIDE
COMPLETE FOR SCREW DRIVER W/O SCREW HOLDING DEVICE
3.5MMO CMF-35000
HANDLE WITH SLEVE CMF-35010
CHEK-HOLDING CMF-35020
TROCAR CMF-35030
DRILL GUIDE CMF-35040
COMPLETE FOR SCREW DRIVER WITH SCREW HOLDING DEVICE
6.0MMO CMF-60000
HANDLE WITH SLEVE CMF-60010
CHEEK – HOLDING CMF-60020
TROCAR CMF-60030
DRILLGUIDE CMF-60040

